

TACTICAL BOUGIETM

Endotracheal tube introducer with Safety Markers



5mm x 70cm

FOR USE WITH ENDOTRACHEAL TUBES 6.0MM ID AND LARGER IN ADULT PATIENTS

COLORED SAFETY MARKINGS | LATEX FREE | SINGLE USE



The Patented safety markings are a unique feature of the disposable Voir Bougie[®]

Adroit Surgical's introducers are manufactured from low density polyethylene which is both firm and flexible for easy insertion.

The Voir Bougie[®] has brightly marked colored bands at both ends that allows rapid accurate insertion to the correct depth by the user and helps to keep the bougie in the same position during use to avoid patient injury or inadvertent removal of the bougie.

Designed specifically for use by EMS professionals.
Fits in paramedic bags, pockets and compact spaces.



Adroit Surgical LLC -13904 Quailbrook Drive Oklahoma City - OK 73134

1-877-7ADROIT
www.adroitsurgical.com

